



Building Leadership Excellence



A Novel Ultrasonic Dynamic Drainage Tester

Jinlan Ju, Xujun Hua and Huining Xiao*



FPIInnovations, Canada

* University of New Brunswick, Canada


May 1-4
PaperCon 2011
Northern Kentucky Convention Center

RETHINK PAPER:
Lean and Green

Water Removal in Papermaking



Cost of Water Removal



Relative cost to remove a given quantity of water:

Former

\$1

Presses

\$58

Dryers

\$220



TAPPI

PaperCon 2011

Conventional Drainage Tester



**Canadian Standard
Freeness Tester:**


**Excellent tool for pulp
quality control**



TAPPI

PaperCon 2011

Dynamic Drainage Jar



A very useful device for determining first pass retention.


Modified Dynamic Drainage Jar




**Dynamic Drainage
Analyzer:**

**A simple test,
but no turbulence during
drainage**

Effect of Mixing on Drainage in Modified DDJ



Fines Distribution of a Sheet Formed in a DDJ Under Continuous Mixing



Other Drainage Testers


Moving Belt Drainage Tester
Pulsed Drainage Device

...




PaperCon 2011


Ultrasonic Dynamic Drainage Tester




Data Collection




Ultrasonic Filtration Screen



Screen Unblocking and Deflocculation by Ultrasound in Standard Retention Test



Drainage Profiles of a Newsprint Furnish At 400 rpm with Mat Formation



TAPPI

PaperCon 2011


Drainage and Retention with Chitosan

At 400 rpm with Mat Formation

Sample	No Mixing	Mixing at 400 rpm	Ultrasound
	Drainage Time (s)		
Pulp Alone	14.1	17.0	11.5
Pulp + Chitosan	12.3	14.1	11.8
	Fines Retention (%)		
Pulp Alone	77.6	79.7	64.1
Pulp + Chitosan	88.6	87.3	80.9
Increase in Retention	11.0	7.6	16.8




Drainage Profiles of a Newsprint Furnish At 1000 rpm with Mat Formation




TAPPI

PaperCon 2011


Retention Is Changed by Ultrasound Intensity




Pressure Profile under Various Turbulence Conditions




Drainage vs. Pressure and Retention




Fines Distribution of a Sheet Formed with Ultrasound




Drainage Mechanism with Ultrasound



Before starting drainage test



Dewatering under mixing



Dewatering under ultrasound

Conclusions

- **The turbulence induced by impellor in a conventional DDJ increases filtration resistance enormously.**
- **Ultrasound can**
 - ✓ **Deflocculate fiber suspension**
 - ✓ **Unblocking forming wire**
 - ✓ **Enlarge dewatering channel in the web**
- **The Ultrasonic DDJ could be very useful laboratory tool to determine drainage under dynamic conditions**
 - ✓ **Drainage profile**
 - ✓ **Vacuum profile**
 - ✓ **Retention**



Acknowledgements

- **T. Owston**
 - ✓ Initial experimental work
- **C. Castro**
 - ✓ Data processing program
- **G. Dorris, G. Laivins, T. Owston & M. Drainville**
 - ✓ Filtrate weight and pressure monitoring
- **I. Tihanyi & N.Vignola**
 - ✓ Vessel and accessory manufacturing
- **M. Laleg, J. Jong, I. Vadeiko, T. Owston and S. Middleton**
 - ✓ Valuable discussion
- **NSERC IPS program**

